









1. Before you start

 Check the contents: MINIPIX device, memory stick with the software, drivers and documentation, Quick Start Guide.



 Warning: the sensor is very fragile – DO NOT touch it or clean it under any circumstances and avoid any mechanical shock!

2. Installation of the PIXETPRO software

- Connect the device to the computer.
- Connect the attached flash to the computer.
- Install the device drivers. On Microsoft Windows run **Drivers/setup.exe** and follow the instructions. On Mac OS X run **Drivers/D2xxHelper_v2.0.0.pkg.** On Linux, it is not necessary install any special drivers.
- Run Pixet/Windows/setup.exe for installation of PIXETPRO on Windows.
- Run the PIXET PRO software from Start Menu shortcut on Windows.
- For information how to use PIXET_{PRO} refer to the corresponding manual in Documentation/Pixet Pro Users Manual.pdf

For questions and support contact us at support@advacam.com

ADVACAM s.r.o.

U Pergamenky 1145/12, CZ 170 00 Praha 7

Czech Republic Tel: +420-603-444112

Email: info@adavcam.com

www.advacam.com

We are environment friendly – all the manuals are available in electronic form on the flash disk.

